

# *BIBO* LED-BITE BLOCK

An intra oral light system which simultaneously, relaxes the patient's jaws as a mouth props and illuminates the working area in the oral cavity.

*Bibo* is extremely simple to use with rechargeable power for the LED light of about 50,000 life hours given light power of 22,000 LUX. Discharge time approximately one hour and charging time about two hours.

An expensive autoclavable bite blocks.


*Bibo* set comes with:

Two small, two large autoclavable silicon bite blocks

One lighting unit with the LED light

One charger

*None illuminated working area*



*Illumination of the working area*



*Wall plug charger*

*Bibo Introductory Price:  
Only \$229.00*

**DENT CORP.**  
Research & Development.

656 Mamaroneck Ave. White Plains, NY 10605  
(800)454-9244/914-682-6600 Fax 914-948-1711  
[www.dentcorp.com](http://www.dentcorp.com) [dentcorp@cloud9.net](mailto:dentcorp@cloud9.net)